
Attenuation of Radiation through the absorption of Gamma Rays



Swapnil Bhatta
The University of Southern Mississippi
March 10th, 2020



Summary of topics:


- Introduction
 - Radiation
 - Gamma Rays
 - History
- Theory
 - Gamma Rays interaction with matter
 - Linear and Mass Attenuation Coefficients
- Experiment
 - Lead
 - Cs-137
 - Geiger-Müller counter
 - Equipment summary
 - Procedure
- Results
- Error and Uncertainty
- Conclusions
- Questions.





Introduction: Radiation


- Emission of energy in the form of EM radiation or particles.
- Occurs naturally in unstable atoms, which radiate characteristic energies to form new stable atoms.
 - Electromagnetic Radiation (microwave, x-rays, gamma rays (γ).
 - Particle Radiation (alpha (α), beta (β), neutrons)
- Ionizing (< 10 electron volts) and nonionizing.






Introduction: History









Wilhelm Röntgen

Henri Becquerel

Ernest Rutherford

Perrie & Marie Curie

First discovery: William Herschel around 1800 (infrared radiation).

Becquerel discovered radioactivity. Marie Curie coined the word, discovered it is specific to given elements of isotopes.

[2]



Introduction: Gamma Rays


- Gamma rays are photons of a very high-frequency electromagnetic wave.
- Generated through unstable atoms radiating energy, or after particle emission
- Nuclear isomer undergoing isomeric transition and releasing photons.
- Discovered by French chemist/physicist Paul Villard.
 - In the early 20th century, detected powerful radiation from radium, distinct from α radiation (Becquerel) and β radiation (Rutherford)






Theory: Gamma Rays and interaction with matter

- Photoelectric effect: Absorbed upon contact, electron ejected.
- Compton Effect: For 100 keV 10 MeV gamma rays, imparts partial energy.
- Pair production: For energies
 >1.02 MeV, gamma ray photons can spontaneously transform into electrons & positrons.





- Given their high energy, gamma rays are quite penetrative. A lot higher than other radiations. Some interact as previously mentioned, some pass through.
- Attenuation of gamma rays is simply the reduction of their intensities.
- Linear attenuation coefficient is a characterization of how easily a material can be penetrated by particles or energy.






The attenuation of gamma rays itself is exponential in nature and follows the equation;

$$I=I_{o}e^{-\mu X}$$

Where I_o is the original intensity of the incoming beam or particles or photons, X is the thickness, and μ is the linear attenuation coefficient.







If half of the incoming intensity is attenuated;

$$\frac{1}{2} I = I_0 e^{-\mu X_{1/2}}$$
or' $\mu = \ln(2)/X_{1/2}$

This gives us two ways to determine the linear attenuation coefficient.









Relationship to various parameters of the absorbing material:

- a) Atomic Number
- b) Density
- c) Thickness
- d) Gamma-Ray energy







Experiment: Lead

Lead is a chemical element with Z=82 denoted as Pb.

- Interesting as it's a cheap and malleable metal.
- Density of high density of 11.34g/cm3.
- Mass attenuation coefficient of lead is expected to be 0.1136 g/cm2 for a 662keV source. [6]






Experiment: Caesium

The gamma ray source: Caesium-137 (Z=55)


- The Caesium nucleus can decay into a more stable nucleus through two different means as shown in the energy level diagram.
- Isomeric transition to a stable Ba-137, by emitting 662 keV of energy.





Experiment: Geiger-Müller counter

- A Geiger-Müller counter is an instrument that detects ionizing radiation.
 - Airtight compartment.
 - Thin conducting wire.
 - High potential difference
 - "avalanche" detection.
 - Needs to be calibrated.
 - Optimal operating voltage.
 - Resolving time.
 - Efficiency.





Experiment: Equipment

- SPECTECH radiation counter
 - High voltage
 - Preset times
 - Preset step-voltages
- Geiger Mueller Counter
- Various radiation sources.
- Lead Absorber Set
 - Lead (0.038-0.25 inches)








Experiment: Procedu

- GM Counter plateau
 - Knee-Voltage
- Background Radiation
- Resolving time
- Efficiency
- Procedures
 - First method
 - Second method
 - Third method
 - Shelfs of the well of the GM tube


Geiger Plataeu




A plateau graph for a Geiger-Müller counter.



Results: Optimal Voltage Issue








Results: Summary of Paramet


- GM Counter plateau
 - Optimal Voltage: 1200V.
- Background Radiation
 - ~30 counts per minute.
- Resolving time
 - 20 μ.s (1-100 μ.s)
- Efficiency
 - 1.33% (0-10%). 0.12% for gamma source.





Results: Procedure One







$$\mu_{1\text{slope}} = 1.09 \pm 0.26 \text{ cm}^{-1}$$
 $X_{1/2}$ (1 (μ/ρ) of 0.096 \pm 0.023 g/cm²

$$\mu_{1\text{slope}} = 1.09 \pm 0.26 \text{ cm}^{-1}$$
 $X_{1/2}$ (1196 counts) was 0.62 cm, $\mu_{1(X1/2)} = 1.12 \text{ cm}^{-1}$ (μ/ρ) of 0.096 \pm 0.023 g/cm² (μ/ρ) of 0.099 g/cm².



Results: Procedure Two







 (μ/ρ) of 0.120 ± 0.015 g/cm² (μ/ρ) of 0.135 g/cm².

 $\mu_{2\text{slope}} = 1.36 \pm 0.17 \text{ cm}^{-1}$ $X_{1/2}$ (452 counts) was 0.453 cm, $\mu_{1(X1/2)} = 1.53 \text{ cm}^{-1}$



Results: Procedure Three





$$\mu_{3\text{slope}} = 1.24 \pm 0.079 \text{ cm}^{-1}$$
 $X_{1/2}$ (452 counts) was 0.429 cm, $\mu_{1(X1/2)} = 1.62 \text{ cm}^{-1}$ (μ/ρ) of 0.109 \pm 0.007 g/cm² (μ/ρ) of 0.142 g/cm².



Results: Summary of Attenuation

- GM Counter plateau
 - Optimal Voltage: 1200V.
- Background Radiation
 - ~30 counts per minute.
- Resolving time
 - 20 μ.s (1-100 μ.s)
- Efficiency
 - 1.33% (0-10%). 0.12% for gamma source.

- Average of μ/ρ (μ from slope).
 - $0.108 \pm 0.029 \text{ g/cm}^2$
- Average of μ/ρ (μ function of half intensity distance).
 - 0.125 g/cm².



Error and Uncertainty

• Using the value of (μ/ρ) obtained from the slopes of $ln(l_0)$ vs Thickness (0.108 ± 0.029 g/cm²) and the accepted mass attenuation coefficient (μ/ρ) of 0.1136 g/cm².

$$\%Error = \frac{|Accepted\ Value - Obtained\ Value|}{Accepted\ Value} \times 100\%$$

$$= \frac{|0.1136 \text{ g/cm2} - 0.108 \text{ g/cm2}|}{0.1136 \text{ g/cm2}} \times 100\%$$

$$=4.93\%$$

 Uncertainty terms in the averaged values through error propagation gives;

$$\Delta(\mu/\rho) = 0.029 \text{ g/cm}2$$

 Using the uncertainty to calculate the percentage fractional uncertainty gives;


$$(4/(\mu p)p)/(400\%) = 26.85\%$$

Withirdh compared to the % error suggests the results are valid.



Conclusion

- Successfully determined the the mass attenuation coefficient of lead through two different approaches.
- Lead is indeed an efficient material to absorb, at least for low energy emissions.
- Could've used more varying density's of Pb as the most accurate procedure did use combinations.
 - The distance between the source and the varying leads could've been as issue.
- Better means of calculating the resolving time as adjusted times had 10-20 counts.
- Better analysis for error propagation.
- Discuss Slide 10 answers.




Citations

- [1] Encyclopedia Britannica. Electromagnetic Spectrum.
- [2] Wikipedia. Paul Villard. https
 https
 https://es.wikipedia.org/wiki/Paul Ulrich Villard#/media/File:Paul Villard.jpg
- [3] WikiCommons. Rutherford.
- [4]https://en.wikipedia.org/wiki/Photoelectric_effect #/media/File:Photoelectric_effect.svg
- [5] Wikibooks. Basic Physics of Nuclear Medicine
- [6] NIST. X-Ray Mass Attenuation Coefficients
- [7] Lead. https://en.wikipedia.org/wiki/Lead
- [11] https://indico.cern.ch/event/593812/contributions/2499768/attachments/14663 50/2268202/KDK Cap Congress 2017 Matthew Stukel.pdf



Thank you.

Questions?



- ▲ Do not perform any work on items labeled with radioactive cautionary signs.
- ⚠ Do not handle items that are labeled with cautionary signs.
- A Do not handle or empty radioactive waste containers.
- ⚠ Do not dust or clean radioactive material work stations.
- A Do not eat, drink or smoke in radioactive material labs.
- ⚠ Do not attemp to clean up a radioactive spill.
- ▲ Wear protective clothing-lab coats, aprons or coveralls.

